

3-Ply Community Face Mask

Product Technical Sheet



DETMOLD
Medical

The **3-Ply Community Mask** is a single-use face mask designed for public or community spaces where social distancing cannot be achieved.

Featuring three layers of lightweight material with two ear loops to keep the mask in place, this mask provides comfort and breathability for the wearer.

The 3-Ply Community Mask can be used in retail, food service and transport environments, as well as for events, festivals and conferences.

This item can be purchased in a carton of 24 boxes of 50 masks, totaling 1,200 masks, or a Shelf Ready Tray containing 30 individual retail packs containing 10 masks in each pack, totaling 300 masks.

This mask **is not** intended for use in clinical environments.




- + Australian Made and Owned
- + Single-use face masks with ear loops
- + Inner comfort layer for comfortable fit
- + Non-sterile and Latex-free

Material Specifications

	Outer Layer	Middle Layer	Inner Layer	Nose Wire
Material Description	Spunbond PP	Meltblown PP	Spunbond PP	Steel
Special Treatment	Hydrophobic	3% Electret	Hydrophilic	Plastic coated
Material gsm	Typically 25gsm	Typically 25-27gsm	Typically 25gsm	3mm wide

Outer Layer

The external hydrophobic layer provides protection for the filtration layer




Middle Layer

The meltblown filtration layer protects from the passage of non-oil based particulates outside the mask.



Inner Layer


Soft and comfortable, the inner hydrophilic (absorbent) layer captures droplets from wearer.



● Droplets from wearer ● Non-oil based particulates

3-Ply Community Face Mask

Product Specification Overview



Product Details

Pack Size	Description	Size	Width	Height	Ear Loop Length	Masks Per Box	Masks Per Carton
10 Pack & SRT	3-Ply Community Face Mask	One Size	175mm ± 2mm	95mm ± 2mm	165mm ± 5mm	10	300
50 Pack	3-Ply Community Face Mask	One Size	175mm ± 2mm	95mm ± 2mm	165mm ± 5mm	50	1200

Instructions for Use

Always check mask fit to ensure it is secured at the bridge of the nose and mouth. There should be minimal gaps between the mask and face. Appropriate hand hygiene should always be used when touching the mask. Masks should not be reapplied and should be discarded safely after use.



1. Hold the mask on the short edges with the ear loop connection points facing outwards.



2. Pinching firmly on the nose piece, place one ear loop around the back of your ear and ensure it is secure.



3. Place the second ear loop up over your other ear and check the mask is evenly aligned on your face.



4. Expand the pleats to ensure a secure and comfortable fit.

5. Check mask fit to ensure it is secured at the bridge of the nose and mouth. There should be minimal gaps between the mask and face.

Storage Recommendations

Store mask in the original packaging, away from contaminated areas, dust, sunlight, extreme temperatures, excessive moisture, and damaging chemicals. When stored in accordance with recommendations, the product is designed for use within 2 years from date of manufacture.



DO NOT USE IF PACKAGING IS DAMAGED



KEEP DRY



DO NOT RE-USE



Caution: The mask is designed for occupational use by healthcare professionals who are properly trained in mask use and limitations. Misuse may result in sickness or death. DO NOT use in any manner not indicated in the Instructions for Use. Always read the label. Discard the mask if it becomes damp, soiled, damaged, or difficult to breathe through. Do not use if the packaging is damaged.

Contact Detmold Medical to place an order today



connect@detmoldmedical.com

1800-DET-MED (1800 338 633)

detmoldmedical.com