

P2 D95 Respirator – FitClip™

The **P2 D95** is a single-use respirator designed to filter out non-oil based airborne particles.

Featuring a lightweight flat fold design, the respirator expands into a convex-shaped device with two bands and a FitClip™ to hold the respirator in place. The knitted head bands are skin friendly and well tolerated.

- + Designed and tested to meet performance requirements of a Class P2 Respirator specified by AS/NZS 1716:2012 Respiratory Protective Devices
- + Non-sterile & Latex Free




Testing Standards

We operate to Quality Management System ISO: 13845. The P2 D95 Respirator is designed and tested to AS/NZS 1716:2012 and offers Class 2 protection against mechanically and thermally generated (dust, mists and fumes) particulates. Raw Materials and Finished Goods are tested routinely in Detmold Medical's validated state-of-the-art laboratory.

All Detmold Medical Respirators are designed to meet Australian conditions, comply with International specifications, and achieve high performance for filtration of sub-micron particles exceeding the 94% and 95% requirements of the AS, NIOSH, and ISO standards, whilst

maintaining a benchmark **low resistance** for breathing. Differential Pressure is the breathability and comfort indicator, which is measured whilst the respirator is worn by determining the difference of pressure across the mask under specific conditions of airflow, temperature, and humidity.




P2 D95 RESPIRATOR - FITCLIP™ PERFORMANCE REQUIREMENTS				
Characteristics	42 CFR Part 84	AS 4381:2015	AS/NZS 1716:2012	Quality
Medical Device Certified System				ISO 13485
ARTG				342752
Particulate Filtration Efficiency (PFE)			Class P2 (>94%)	
Bacterial Filtration Efficiency (BFE)				
Blood Splash Resistance				

Note: These performance requirements are measured on unworn respirators. The respirator may not meet these performance requirements

during wear due to individual fit characteristics. Please consult your fit testing procedure. This product is designed for use in a healthcare setting.

Material Specifications



P2 D95 Respirator - FitClip™




Product Details

Item No.	Description	Size	Qty Per Box	Qty Per Carton	Shipper Weight	Shipper dimensions
W887M	P2 D95 Respirator with regular size FitClip™	X-Small	25	600	5.5kg	324 x 505 x 333mm
W575M		Small				
W376M		Regular				
W486M		Large				
W187M		X-Large				

Item No.	Description	Size	Inner Quantity	Carton Quantity	Shipper Weight	Shipper dimensions
Supplied with all sizes	FitClip™	Regular (35 x 20mm)	N/A	N/A	N/A	N/A
W008S9001		X-Large (69 x 19mm)	50 pieces	N/A	N/A	N/A
W008S0001		X-Large (69 x 19mm)	N/A	1000 pieces	1.2kg	125 x 176 x 150mm

Securing FitClip™



Securing the Class P2 D95 Respirator with FitClip™

The P2 D95 Respirator is designed to achieve a comfortable and adequate seal by joining the respirator's two bands with a FitClip™. The Detmold Medical Respirators have the Regular FitClip™ already attached, and where a longer band is required, the XL FitClip™ can be ordered separately. An inadequate seal between the Respirator and the face can result in contaminated air entering the respirator. Please refer to your appropriate Respirator Fit Test Procedure to determine the required combination of respirator and FitClip™.

Instructions for Use

Applying the Respirator | Prior to donning or doffing, ensure hand hygiene is undertaken.

There should be no gaps between the face and respirator, providing an adequate seal[^]. Once the respirator is applied it should not be touched or reapplied. **Please discard safely after use.**



1. Move the FitClip™ to the band on your preferred side for donning*.
2. Open the mask, flattening the nose piece.



3. Lean chin into the mask and lift bands toward the crown of the head securing onto the FitClip™.
4. Adjust mask bands for a secure and comfortable fit.



5. Using both hands, firmly mould nose piece over bridge of the nose.
6. Pinch the nose piece (angled over the nostril) and drop the jaw gently, allowing the mask to fit to the natural shape of your face.



7. Look left to right, then up and down ensuring there is no movement in the mask at the nose piece.
8. Inhale and exhale deeply checking for escape of air. If air is escaping, refit the mask.

Removing the Mask



1. Without touching the mask, follow bands towards FitClip™ and gently unhook, keeping tension on the bands.
2. Lean forward maintaining control of the band tension, whilst slowly removing the mask away from the face. Ensure the mask does not flick forward or back onto the face.
3. Discard responsibly.

Storage Recommendations

Store respirator in the original packaging, away from contaminated areas, dust, sunlight, extreme temperatures, excessive moisture, and damaging chemicals. Storage temperature -20°C to +30°C at <80% relative humidity. When stored in accordance with recommendations, the product is designed for use within 5 years (2 years for W887M), from date of manufacture.



DO NOT USE IF PACKAGING IS DAMAGED



KEEP DRY



DO NOT RE-USE



Caution: The respirator is designed for occupational use by healthcare professionals who are properly trained in respirator use and limitations. Misuse may result in sickness or death. DO NOT use in any manner not indicated in the Instructions for Use. Always read the label. Discard the respirator if it becomes damp, soiled, damaged, or difficult to breathe through. Do not use if the packaging is damaged.

Contact Detmold Medical to place an order today

©2024 Detmold Medical. All Rights Reserved. February 2024. #DM-RP2-2024-002

*FitClip™ design may vary. Contact Detmold Medical if an alternate FitClip™ size is required. *Follow your relevant fit check test procedure.



connect@detmoldmedical.com

1800 338 633

detmoldmedical.com