

RETAIL PACK | 20 PIECES
EXTRA SMALL | HEAD BANDS
N95 SURGICAL RESPIRATOR
MASK

The Detmold Medical EXTRA SMALL N95 Surgical Respirator Mask is a single-use surgical respirator to protect those working in hospitals, healthcare, and clinical environments where a person is at risk of blood or fluid splash.

The N95 mask filters out non-oil-based airborne particulates and bacteria. Please see the adjacent Features section for this N95 Surgical Respirator Mask's performance.

What's In The Retail Mask ...

[More info online](#)



Item Code	W125M9001
Range	N95 Surgical Respirator Masks
Colour	White
Carton	20
Pieces Per Sleeve	20

Carton Width (mm)	125
Carton Length (mm)	165
Carton Height (mm)	105
Carton Gross Weight (kg)	0.2
MOQ - Stock	20

Manufacturing certifications	This product is made in an ISO: 13485 medical device manufacturing facility.
------------------------------	--