

Respiratory PPE
Selection
Guide



Page Number	13	12	15	18	19	20	21	16
Protection		L3BFEARTG	L3N95BFEARTG	L3N95BFEARTG	N95P2BFEARTG	P2ARTG	P2	P2ARTG
Sizes	One Size	One Size	XS - XL	One Size	One Size	One Size	One Size	XS - XL
Fastener	Standard Ear Loops	Standard Ear Loops / Short Ear Loops / Head Tie / Visor	Headband / FitClip™	Headband	Headband	Headband	Headband	Headband / FitClip™
Healthcare		Operating Rooms (OR) Surgeries Emergency Rooms ER Frontline Health incl. ambulance, ambulatory care Haematology Oncology Infection Control Patient Care Hospitals GPs	Operating Rooms (OR) Surgeries Emergency Rooms ER Frontline Health incl. ambulance, ambulatory care Haematology Oncology Infection Control Patient Care Hospitals GPs	Operating Rooms (OR) Surgeries Emergency Rooms ER Frontline Health incl. ambulance, ambulatory care Haematology Oncology Infection Control Patient Care Hospitals GPs	Infection Control Patient Care Hospitals GPs			
Veterinary			Vet Surgeries / Animal Care	Vet Surgeries / Animal Care	Vet Surgeries / Animal Care			
Aged Care				Aged Care				
Beauty Therapy	Beauty			Beauty	Beauty			Beauty
Crowded Environments	Events / Travel Public / Retail / Community					Events / Travel Public / Retail / Community		Events / Travel Public / Retail / Community
Industrial						Construction / Demolition Farming (Grain / Poultry) Mining / O&G Manufacturing Industrial	Construction / Demolition Farming (Grain / Poultry) Mining / O&G Manufacturing Industrial	Construction / Demolition Farming (Grain / Poultry) Mining / O&G Manufacturing Industrial

Protection Levels:

- N95N95 - Particulate Filtration Efficiency (NIOSH 42CFR, Part 84)
- P2Class P2 - Particulate Filtration Efficiency (AS/NZS 1716:2012)

L3L3 - Blood/Fluid Resistance (AS4381:2015)

BFEBacterial Filtration Efficiency (AS4381:2015)

ARTGIncluded on the Australian Register of Therapeutic Goods (ARTG)

